

For Research Use Only

GSTT2B Recombinant antibody, PBS Only (Detector)



Catalog Number: 83469-1-PBS

Basic Information

Catalog Number:

83469-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11828

GenBank Accession Number:

BC071700

GeneID (NCBI):

653689

UNIPROT ID:

POCG30

Full Name:

glutathione S-transferase theta 2B
(gene/pseudogene)

Calculated MW:

244 aa, 28 kDa

Purification Method:

Protein A purification

CloneNo.:

240493D2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00473-1, GSTT2B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83469-2-PBS. Detection antibody: 83469-1-PBS. Standard: Ag11828. Range: 0.78-100 ng/mL