

For Research Use Only

GSTT2B Recombinant antibody, PBS Only (Capture)



Catalog Number: 83469-2-PBS

Basic Information

Catalog Number: 83469-2-PBS	GenBank Accession Number: BC071700	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 653689	CloneNo.: 240493E5
Source: Rabbit	UNIPROT ID: POCG30	
Isotype: IgG	Full Name: glutathione S-transferase theta 2B (gene/pseudogene)	
Immunogen Catalog Number: AG11828	Calculated MW: 244 aa, 28 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00473-1, GSTT2B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83469-2-PBS. Detection antibody: 83469-1-PBS. Standard: Ag11828. Range: 0.78-100 ng/mL