

For Research Use Only

FAM116A Recombinant antibody, PBS Only (Capture)



Catalog Number: 83473-1-PBS

Basic Information

Catalog Number: 83473-1-PBS	GenBank Accession Number: BC040291	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 201627	CloneNo.: 240486A7
Source: Rabbit	UNIPROT ID: Q8IWF6	
Isotype: IgG	Full Name: family with sequence similarity 116, member A	
Immunogen Catalog Number: AG22052	Calculated MW: 608 aa, 70 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00477-1, FAM116A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83473-1-PBS. Detection antibody: 83473-3-PBS. Standard: Ag22052. Range: 0.78-100 ng/mL.



Cytometric bead array standard curve of MP00477-2, FAM116A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83473-1-PBS. Detection antibody: 83473-2-PBS. Standard: Ag22052. Range: 0.78-100 ng/mL.