

For Research Use Only

MIG6; ERRFI1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83475-3-PBS

Basic Information

Catalog Number: 83475-3-PBS	GenBank Accession Number: NM_018948	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 54206	CloneNo.: 240447A4
Source: Rabbit	UNIPROT ID: Q9UJM3	
Isotype: IgG	Full Name: ERBB receptor feedback inhibitor 1	
Immunogen Catalog Number: AG33341	Calculated MW: 51kd	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00464-2, MIG6; ERRF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83475-1-PBS. Detection antibody: 83475-3-PBS. Standard: Ag33341. Range: 0.625-80 ng/mL.