

For Research Use Only

FOXC2 Recombinant antibody, PBS Only (Capture)

Catalog Number:83476-7-PBS



Basic Information

Catalog Number:

83476-7-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33121

GenBank Accession Number:

NM_005251

GeneID (NCBI):

2303

UNIPROT ID:

Q99958

Full Name:

forkhead box C2 (MFH-1,
mesenchyme forkhead 1)

Calculated MW:

54kd

Purification Method:

Protein A purification

CloneNo.:

240403B7

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00470-4, Human FOXC2 Recombinant Matched Antibody Pair - PBS only. 83476-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag33121. 83476-5-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean FOXC2 concentration was determined to be 2.8 ng/mL in HeLa cell extract based on a 1.5 mg/mL extract load.