

For Research Use Only

C2orf7 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83478-1-PBS

Basic Information

Catalog Number: 83478-1-PBS	GenBank Accession Number: BC005069	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 84279	CloneNo.: 240494A9
Source: Rabbit	UNIPROT ID: Q9B5G0	
Isotype: IgG	Full Name: chromosome 2 open reading frame 7	
Immunogen Catalog Number: AG14678	Calculated MW: 188 aa, 21 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00468-2, C2orf7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83478-3-PBS. Detection antibody: 83478-1-PBS. Standard: Ag14678. Range: 1.25-80 ng/mL