

For Research Use Only

TRIM22 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83479-1-PBS

Basic Information

Catalog Number:

83479-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4683

GenBank Accession Number:

BC035582

GeneID (NCBI):

10346

UNIPROT ID:

Q8IYM9

Full Name:

tripartite motif-containing 22

Calculated MW:

498 aa, 57 kDa

Purification Method:

Protein A purification

CloneNo.:

240411A12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00462-1, TRIM22 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83479-1-PBS. Detection antibody: 83479-3-PBS. Standard: Ag4683. Range: 0.625-80 ng/mL