

For Research Use Only

# HMGB2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83482-2-PBS

## Basic Information

Catalog Number: 83482-2-PBS	GenBank Accession Number: BC001063	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3148	CloneNo.: 240453G9
Source: Rabbit	UNIPROT ID: P26583	
Isotype: IgG	Full Name: high-mobility group box 2	
Immunogen Catalog Number: AG6135	Calculated MW: 24 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00480-1, HMGB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83482-2-PBS. Detection antibody: 83482-1-PBS. Standard: Ag6135. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP00480-2, HMGB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83482-2-PBS. Detection antibody: 83482-3-PBS. Standard: Ag6135. Range: 3.125-100 ng/mL.