

For Research Use Only

HMGB2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83482-3-PBS

Basic Information

Catalog Number: 83482-3-PBS	GenBank Accession Number: BC001063	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3148	CloneNo.: 240453F8
Source: Rabbit	UNIPROT ID: P26583	
Isotype: IgG	Full Name: high-mobility group box 2	
Immunogen Catalog Number: AG6135	Calculated MW: 24 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00480-2, HMGB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83482-2-PBS. Detection antibody: 83482-3-PBS. Standard: Ag6135. Range: 3.125-100 ng/mL.