

For Research Use Only

SYNJ2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83483-1-PBS



Basic Information

| | | |
|--|--|---|
| Catalog Number: 83483-1-PBS | GenBank Accession Number: BC043277 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 8871 | CloneNo.: 240456E11 |
| Source: Rabbit | UNIPROT ID: O15056 | |
| Isotype: IgG | Full Name: synaptojanin 2 | |
| Immunogen Catalog Number: AG4855 | Calculated MW: 1496 aa, 166 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00479-1, SYNJ2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83483-1-PBS. Detection antibody: 83483-3-PBS. Standard: Ag4855. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00479-2, SYNJ2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83483-1-PBS. Detection antibody: 83483-2-PBS. Standard: Ag4855. Range: 0.313-40 ng/mL



Bi-layer interferometry (BLI) kinetic assay of 83483-1-PBS against Human SYNJ2 was performed. The affinity constant is 0.706 nM.