

For Research Use Only

SYNJ2 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83483-2-PBS



Basic Information

Catalog Number: 83483-2-PBS	GenBank Accession Number: BC043277	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8871	CloneNo.: 240456G1
Source: Rabbit	UNIPROT ID: O15056	
Isotype: IgG	Full Name: synaptojanin 2	
Immunogen Catalog Number: AG4855	Calculated MW: 1496 aa, 166 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00479-2, SYNJ2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83483-1-PBS. Detection antibody: 83483-2-PBS. Standard: Ag4855. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00479-3, SYNJ2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83483-2-PBS. Detection antibody: 83483-3-PBS. Standard: Ag4855. Range: 1.25-40 ng/mL



Bi-layer interferometry (BLI) kinetic assay of 83483-2-PBS against Human SYNJ2 was performed. The affinity constant is below 1 pM.