

For Research Use Only

Mouse MMR/CD206 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83485-3-PBS

Basic Information

Catalog Number:

83485-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_008625.2

GeneID (NCBI):

17533

UNIPROT ID:

Q2HZ94

Full Name:

mannose receptor, C type 1

Calculated MW:

165 kDa

Purification Method:

Protein A purification

CloneNo.:

240642G1

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00623-1, MOUSE MMR/CD206 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83485-4-PBS. Detection antibody: 83485-3-PBS. Standard: Eg0790. Range: 1.563-100 ng/mL.