

For Research Use Only

FAM38B Recombinant antibody, PBS Only (Capture)



Catalog Number: 83488-4-PBS

Basic Information

Catalog Number: 83488-4-PBS	GenBank Accession Number: AB527139	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 63895	CloneNo.: 240406E5
Source: Rabbit	UNIPROT ID: Q9H5I5	
Isotype: IgG	Full Name: family with sequence similarity 38, member B	
Immunogen Catalog Number: AG24528	Calculated MW: 318 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00488-2, FAM38B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83488-4-PBS. Detection antibody: 83488-1-PBS. Standard: Ag24528. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00488-3, FAM38B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83488-4-PBS. Detection antibody: 83488-2-PBS. Standard: Ag24528. Range: 0.313-40 ng/mL.