

For Research Use Only

ARID4A Recombinant antibody, PBS Only (Detector)

Catalog Number:83489-4-PBS



Basic Information

Catalog Number: 83489-4-PBS	GenBank Accession Number: NM_002892.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5926	CloneNo.: 240436G5
Source: Rabbit	UNIPROT ID: P29374	
Isotype: IgG	Full Name: AT rich interactive domain 4A (RBP1-like)	
Immunogen Catalog Number: AG33114	Calculated MW: 143 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00489-1, ARID4A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83489-1-PBS. Detection antibody: 83489-4-PBS. Standard: Ag33114. Range: 0.313-40 ng/mL