

For Research Use Only

NELF Recombinant antibody, PBS Only (Capture)

Catalog Number: 83490-3-PBS



Basic Information

Catalog Number: 83490-3-PBS	GenBank Accession Number: BC016862	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 26012	CloneNo.: 240482C8
Source: Rabbit	UNIPROT ID: Q6X4W1	
Isotype: IgG	Full Name: nasal embryonic LHRH factor	
Immunogen Catalog Number: AG34902	Calculated MW: 507 aa, 56 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00490-2, NELF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83490-3-PBS. Detection antibody: 83490-2-PBS. Standard: Ag34902. Range: 0.625-80 ng/mL