

For Research Use Only

C6orf130 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83500-1-PBS

Basic Information

Catalog Number: 83500-1-PBS	GenBank Accession Number: BC011709	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 221443	CloneNo.: 240455A2
Source: Rabbit	UNIPROT ID: Q9Y530	
Isotype: IgG	Full Name: chromosome 6 open reading frame 130	
Immunogen Catalog Number: AG19341	Calculated MW: 152 aa, 17 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00483-1, C6orf130 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83500-2-PBS. Detection antibody: 83500-1-PBS. Standard: Ag19341. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP00483-3, C6orf130 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83500-1-PBS. Detection antibody: 83500-3-PBS. Standard: Ag19341. Range: 6.25-100 ng/mL