

For Research Use Only

# HEY2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83501-2-PBS

## Basic Information

<b>Catalog Number:</b> 83501-2-PBS	<b>GenBank Accession Number:</b> BC007707	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 23493	<b>CloneNo.:</b> 240474B2
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9UBP5	
<b>Isotype:</b> IgG	<b>Full Name:</b> hairy/enhancer-of-split related with YRPW motif 2	
<b>Immunogen Catalog Number:</b> AG30988	<b>Calculated MW:</b> 36 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00482-2, HEY2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83501-2-PBS. Detection antibody: 83501-4-PBS. Standard: Ag30988. Range: 0.625-40 ng/mL