

For Research Use Only

# HEY2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83501-3-PBS

## Basic Information

Catalog Number:

83501-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30988

GenBank Accession Number:

BC007707

GeneID (NCBI):

23493

UNIPROT ID:

Q9UBP5

Full Name:

hairy/enhancer-of-split related with YRPW motif 2

Calculated MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

240474C3

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00482-1, HEY2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83501-3-PBS. Detection antibody: 83501-5-PBS. Standard: Ag30988. Range: 0.625-40 ng/mL