

For Research Use Only

APOA1 Recombinant antibody, PBS Only



Catalog Number: 83504-1-PBS

Basic Information

Catalog Number:

83504-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC005380

GeneID (NCBI):

335

UNIPROT ID:

P02647

Full Name:

apolipoprotein A-I

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

240269F1

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 83504-1-RR (Apolipoprotein A1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83504-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83504-1-RR against Human Apolipoprotein A1 were performed. The affinity constant is 4.05 nM.