

For Research Use Only

C10orf97 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83508-1-PBS

Basic Information

Catalog Number: 83508-1-PBS	GenBank Accession Number: BC067799	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 80013	CloneNo.: 240508A7
Source: Rabbit	UNIPROT ID: Q9H8M7	
Isotype: IgG	Full Name: chromosome 10 open reading frame 97	
Immunogen Catalog Number: AG23520		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00500-1, C10orf97 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83508-1-PBS. Detection antibody: 83508-2-PBS. Standard: Ag23520. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00500-2, C10orf97 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83508-1-PBS. Detection antibody: 83508-3-PBS. Standard: Ag23520. Range: 0.625-40 ng/mL.