

For Research Use Only

C10orf97 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83508-2-PBS

Basic Information

Catalog Number:

83508-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23520

GenBank Accession Number:

BC067799

GeneID (NCBI):

80013

UNIPROT ID:

Q9H8M7

Full Name:

chromosome 10 open reading frame

97

Purification Method:

Protein A purification

CloneNo.:

240508D4

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00500-1, C10orf97 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83508-1-PBS. Detection antibody: 83508-2-PBS. Standard: Ag23520. Range: 0.313-40 ng/mL.