

For Research Use Only

FBXO38 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83509-2-PBS

Basic Information

Catalog Number: 83509-2-PBS	GenBank Accession Number: XM_005268513.1	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 81545	CloneNo.: 240513G2
Source: Rabbit	UNIPROT ID: Q6PIJ6	
Isotype: IgG	Full Name: F-box protein 38	
Immunogen Catalog Number: AG35444	Calculated MW: 134kDa, 1188aa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00502-1, FBXO38 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83509-2-PBS. Detection antibody: 83509-3-PBS. Standard: Ag35444 Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP00502-2, FBXO38 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83509-4-PBS. Detection antibody: 83509-2-PBS. Standard: Ag35444. Range: 0.781-100 ng/mL