

For Research Use Only

FBXO38 Recombinant antibody, PBS Only



Catalog Number: 83509-5-PBS

Basic Information

Catalog Number: 83509-5-PBS	GenBank Accession Number: XM_005268513.1	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 81545	CloneNo.: 240513A4
Source: Rabbit	UNIPROT ID: Q6PIJ6	
Isotype: IgG	Full Name: F-box protein 38	
Immunogen Catalog Number: AG35444	Calculated MW: 134kDa, 1188aa	

Applications

Tested Applications:
ELISA, IF/ICC

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using FBXO38 antibody (83509-5-RR, Clone: 240513A4) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1), CL594-Phalloidin (red). This data was developed using the same antibody clone with 83509-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83509-5-RR against Human FBXO38 were performed. The affinity constant is 0.233 nM.