

For Research Use Only

GPR105 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83511-1-PBS

Basic Information

Catalog Number:

83511-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14105

GenBank Accession Number:

BC034989

GeneID (NCBI):

9934

UNIPROT ID:

Q15391

Full Name:

purinergic receptor P2Y, G-protein coupled, 14

Calculated MW:

338 aa, 39 kDa

Purification Method:

Protein A purification

CloneNo.:

240498B10

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00509-1, GPR105 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83511-2-PBS. Detection antibody: 83511-1-PBS. Standard: Ag14105. Range: 6.25-100 ng/mL