

For Research Use Only

# Mouse IL-3 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83512-1-PBS

## Basic Information

Catalog Number: 83512-1-PBS	GenBank Accession Number: NM_010556.4	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 16187	CloneNo.: 240231A11
Source: Rabbit	UNIPROT ID: P01586	
Isotype: IgG	Full Name: interleukin 3	
	Calculated MW: 18KD	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Mouse

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00510-1, MOUSE IL-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83512-4-PBS. Detection antibody: 83512-1-PBS. Standard: Eg0668. Range: 0.313-20 ng/mL.



Cytometric bead array standard curve of MP00510-2, MOUSE IL-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83512-1-PBS. Detection antibody: 83512-3-PBS. Standard: Eg0668. Range: 0.625-20 ng/mL.