

For Research Use Only

Mouse IL-3 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83512-3-PBS

Basic Information

Catalog Number: 83512-3-PBS	GenBank Accession Number: NM_010556.4	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 16187	CloneNo.: 240231F6
Source: Rabbit	UNIPROT ID: P01586	
Isotype: IgG	Full Name: interleukin 3	
	Calculated MW: 18KD	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00510-2, MOUSE IL-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83512-1-PBS. Detection antibody: 83512-3-PBS. Standard: Eg0668. Range: 0.625-20 ng/mL.



Cytometric bead array standard curve of MP00510-3, MOUSE IL-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83512-3-PBS. Detection antibody: 83512-2-PBS. Standard: Eg0668. Range: 0.625-80 ng/mL.