

For Research Use Only

CD4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83513-1-PBS

Basic Information

Catalog Number: 83513-1-PBS	GenBank Accession Number: BC025782	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 920	CloneNo.: 240427A8
Source: Rabbit	ENSEMBL Gene ID: ENSG00000010610	
Isotype: IgG	UNIPROT ID: P01730	
	Full Name: CD4 molecule	
	Calculated MW: 55 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00503-1, CD4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83513-1-PBS. Detection antibody: 83513-5-PBS. Standard: Eg0333. Range: 0.313-40 ng/mL