

For Research Use Only

MRPL40 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83515-1-PBS

Basic Information

Catalog Number: 83515-1-PBS	GenBank Accession Number: BC009707	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 64976	CloneNo.: 240506A8
Source: Rabbit	UNIPROT ID: Q9NQ50	
Isotype: IgG	Full Name: mitochondrial ribosomal protein L40	
Immunogen Catalog Number: AG10341	Calculated MW: 206 aa, 24 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00514-1, MRPL40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83515-1-PBS. Detection antibody: 83515-3-PBS. Standard: Ag10341. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00514-2, MRPL40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83515-2-PBS. Detection antibody: 83515-1-PBS. Standard: Ag10341. Range: 3.13-100 ng/mL