

For Research Use Only

# Mouse PD-1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83516-1-PBS

## Basic Information

Catalog Number:

83516-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_008798.2

GeneID (NCBI):

18566

UNIPROT ID:

Q02242

Full Name:

programmed cell death 1

Calculated MW:

32kDa

Purification Method:

Protein A purification

CloneNo.:

240390B4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00496-3, MOUSE PD-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83516-3-PBS. Detection antibody: 83516-1-PBS. Standard: Eg0918. Range: 0.156-20 ng/mL