

For Research Use Only

Mouse PD-1 Recombinant antibody, PBS Only (Capture)



Catalog Number:83516-3-PBS

Basic Information

Catalog Number:

83516-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_008798.2

GeneID (NCBI):

18566

UNIPROT ID:

Q02242

Full Name:

programmed cell death 1

Calculated MW:

32kDa

Purification Method:

Protein A purification

CloneNo.:

240389B8

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00496-3, MOUSE PD-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83516-3-PBS. Detection antibody: 83516-1-PBS. Standard: Eg0918. Range: 0.156-20 ng/mL