

For Research Use Only

SOD2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83519-1-PBS

Basic Information

Catalog Number: 83519-1-PBS	GenBank Accession Number: BC016934	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6648	CloneNo.: 240462C5
Source: Rabbit	UNIPROT ID: P04179	
Isotype: IgG	Full Name: superoxide dismutase 2, mitochondrial	
Immunogen Catalog Number: AG21388	Calculated MW: 25 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00517-1, SOD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83519-2-PBS. Detection antibody: 83519-1-PBS. Standard: Ag21388. Range: 1.56-200 ng/mL.