

For Research Use Only

SMPD2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83520-2-PBS

Basic Information

| | | |
|--|---|---|
| Catalog Number: 83520-2-PBS | GenBank Accession Number: BC000038 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 6610 | CloneNo.: 240400D4 |
| Source: Rabbit | UNIPROT ID: O60906 | |
| Isotype: IgG | Full Name: sphingomyelin phosphodiesterase 2, neutral membrane (neutral sphingomyelinase) | |
| Immunogen Catalog Number: AG7367 | Calculated MW: 48 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00518-2, SMPD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83520-2-PBS. Detection antibody: 83520-4-PBS. Standard: Ag7367. Range: 0.78-100 ng/mL.