

For Research Use Only

CD86 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83523-1-PBS

Basic Information

Catalog Number: 83523-1-PBS	GenBank Accession Number: NM_175862	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 942	CloneNo.: 240429B10
Source: Rabbit	ENSEMBL Gene ID: ENSG00000114013	
Isotype: IgG	UNIPROT ID: P42081	
	Full Name: CD86 molecule	
	Calculated MW: 329 aa, 38 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00511-1, CD86 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83523-3-PBS. Detection antibody: 83523-1-PBS. Standard: Eg0966. Range: 0.313-40 ng/mL.