

For Research Use Only

CEACAM1 Recombinant antibody, PBS Only (Detector)



Catalog Number:83524-2-PBS

Basic Information

Catalog Number:

83524-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC014473

GeneID (NCBI):

634

UNIPROT ID:

P13688

Full Name:

carcinoembryonic antigen-related
cell adhesion molecule 1 (biliary
glycoprotein)

Calculated MW:

526 aa, 58 kDa

Purification Method:

Protein A purification

CloneNo.:

240497H8

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00520-3, CEACAM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83524-3-PBS. Detection antibody: 83524-2-PBS. Standard: Eg0314. Range: 0.313-40 ng/mL.