

For Research Use Only

Rat Adipoq Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83525-3-PBS

Basic Information

Catalog Number: 83525-3-PBS	GenBank Accession Number: NM_144744.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 246253	CloneNo.: 240392C9
Source: Rabbit	UNIPROT ID: Q8K3R4	
Isotype: IgG	Full Name: adiponectin, C1Q and collagen domain containing	
	Calculated MW: 26 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Rat

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00522-1, RAT Adiponectin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83525-3-PBS. Detection antibody: 83525-2-PBS. Standard: Eg0872. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00522-3, RAT Adiponectin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83525-1-PBS. Detection antibody: 83525-3-PBS. Standard: Eg0872. Range: 0.625-80 ng/mL.