

For Research Use Only

SLC35B4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83527-1-PBS

Basic Information

Catalog Number:

83527-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22669

GenBank Accession Number:

BC008413

GeneID (NCBI):

84912

UNIPROT ID:

Q96950

Full Name:

solute carrier family 35, member B4

Purification Method:

Protein A purification

CloneNo.:

240332E12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00497-2, SLC35B4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83527-1-PBS. Detection antibody: 83527-2-PBS. Standard: Ag22669. Range: 2.5-80 ng/mL.