

For Research Use Only

TRIM21 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83530-2-PBS

Basic Information

Catalog Number: 83530-2-PBS	GenBank Accession Number: BC010861	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6737	CloneNo.: 240491D2
Source: Rabbit	UNIPROT ID: P19474	
Isotype: IgG	Full Name: tripartite motif-containing 21	
Immunogen Catalog Number: AG28377	Calculated MW: 475 aa, 54 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00505-1, TRIM21 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83530-2-PBS. Detection antibody: 83530-1-PBS. Standard: Ag28377. Range: 1.56-100 ng/mL