

For Research Use Only

METTL9 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83531-3-PBS

Basic Information

| | | |
|-------------------------------------|--|--|
| Catalog Number: 83531-3-PBS | GenBank Accession Number: BC000195 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 51108 | CloneNo.: 240510A10 |
| Source: Rabbit | UNIPROT ID: Q9H1A3 | |
| Isotype: IgG | Full Name: methyltransferase like 9 | |
| Immunogen Catalog Number: AG7191 | Calculated MW: 37 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00508-1, METTL9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83531-4-PBS. Detection antibody: 83531-3-PBS. Standard: Ag7191. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00508-3, METTL9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83531-2-PBS. Detection antibody: 83531-3-PBS. Standard: Ag7191. Range: 0.625-80 ng/mL