

For Research Use Only

GBP3 Recombinant antibody, PBS Only



Catalog Number: 83534-1-PBS

Basic Information

Catalog Number: 83534-1-PBS	GenBank Accession Number: BC140837	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2635	CloneNo.: 230178D8
Source: Rabbit	UNIPROT ID: Q9HOR5	
Isotype: IgG	Full Name: guanylate binding protein 3	
Immunogen Catalog Number: AG33889		

Applications

Tested Applications:
ELISA, WB

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83534-1-RR (Gbp3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83534-1-PBS in a different storage buffer formulation.



1×10^6 HT-29 cells were intracellularly stained with 0.25 ug Gbp3 Recombinant Antibody (83534-1-RR, Clone:230178D8) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)(red), or 0.25 ug Isotype Control (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C). This data was developed using the same antibody clone with 83534-1-PBS in a different storage buffer formulation.