

For Research Use Only

IL-3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83544-2-PBS

Basic Information

Catalog Number: 83544-2-PBS	GenBank Accession Number: NM_000588.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3562	CloneNo.: 240417E1
Source: Rabbit	UNIPROT ID: P08700	
Isotype: IgG	Full Name: interleukin 3 (colony-stimulating factor, multiple)	
	Calculated MW: 17 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00526-2, IL-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83544-3-PBS. Detection antibody: 83544-2-PBS. Standard: Eg0829. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00526-3, IL-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83544-4-PBS. Detection antibody: 83544-2-PBS. Standard: Eg0829. Range: 0.313-20 ng/mL