

For Research Use Only

MBOAT7 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83546-1-PBS

Basic Information

Catalog Number: 83546-1-PBS	GenBank Accession Number: BC006309	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 79143	CloneNo.: 240505B1
Source: Rabbit	UNIPROT ID: Q96N66	
Isotype: IgG	Full Name: membrane bound O-acyltransferase domain containing 7	
Immunogen Catalog Number: AG24466	Calculated MW: 472 aa, 53 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00530-1, MBOAT7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83546-1-PBS. Detection antibody: 83546-2-PBS. Standard: Ag24466. Range: 1.56-100 ng/mL.