

For Research Use Only

Mouse Thrombopoietin Recombinant antibody, PBS Only (Capture)



Catalog Number: 83549-1-PBS

Basic Information

Catalog Number:

83549-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001173505.1

GeneID (NCBI):

21832

UNIPROT ID:

P40226

Full Name:

thrombopoietin

Calculated MW:

38KD

Purification Method:

Protein A purification

CloneNo.:

240203A2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00544-1, MOUSE Thrombopoietin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83549-1-PBS. Detection antibody: 83549-2-PBS. Standard: Eg0678. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00544-3, MOUSE Thrombopoietin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83549-1-PBS. Detection antibody: 83549-4-PBS. Standard: Eg0678. Range: 0.156-20 ng/mL