

For Research Use Only

GPR161 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83557-5-PBS

Basic Information

Catalog Number:

83557-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30998

GenBank Accession Number:

BC028163

GeneID (NCBI):

23432

UNIPROT ID:

Q8N6U8

Full Name:

G protein-coupled receptor 161

Calculated MW:

529 aa, 59 kDa

Purification Method:

Protein A purification

CloneNo.:

240522F12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00542-1, GPR161 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83557-5-PBS. Detection antibody: 83557-2-PBS. Standard: Ag30998. Range: 0.156- 20 ng/mL