

For Research Use Only

S100A14 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83559-4-PBS

Basic Information

Catalog Number: 83559-4-PBS	GenBank Accession Number: BC005019	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57402	CloneNo.: 240528G11
Source: Rabbit	UNIPROT ID: Q9HCY8	
Isotype: IgG	Full Name: S100 calcium binding protein A14	
Immunogen Catalog Number: AG0757	Calculated MW: 12 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00547-2, S100A14 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83559-1-PBS. Detection antibody: 83559-4-PBS. Standard: Ag0757. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP00547-3, S100A14 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83559-4-PBS. Detection antibody: 83559-3-PBS. Standard: Ag0757. Range: 0.781-100 ng/mL