

For Research Use Only

# MAGI2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83561-1-PBS

## Basic Information

Catalog Number:

83561-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18279

GenBank Accession Number:

BC150277

GeneID (NCBI):

9863

UNIPROT ID:

Q86UL8

Full Name:

membrane associated guanylate kinase, WW and PDZ domain containing 2

Calculated MW:

1455 aa, 159 kDa

Purification Method:

Protein A purification

CloneNo.:

240512D2

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00529-3, MAGI2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83561-2-PBS. Detection antibody: 83561-1-PBS. Standard: Ag18279. Range: 0.781-100 ng/mL.