

For Research Use Only

MAGI2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83561-2-PBS

Basic Information

Catalog Number: 83561-2-PBS	GenBank Accession Number: BC150277	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9863	CloneNo.: 240512G6
Source: Rabbit	UNIPROT ID: Q86UL8	
Isotype: IgG	Full Name: membrane associated guanylate kinase, WW and PDZ domain containing 2	
Immunogen Catalog Number: AG18279	Calculated MW: 1455 aa, 159 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00529-2, MAGI2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83561-2-PBS. Detection antibody: 83561-4-PBS. Standard: Ag18279. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00529-3, MAGI2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83561-2-PBS. Detection antibody: 83561-1-PBS. Standard: Ag18279. Range: 0.781-100 ng/mL.