

For Research Use Only

MAGI2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83561-3-PBS

Basic Information

Catalog Number:

83561-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18279

GenBank Accession Number:

BC150277

GeneID (NCBI):

9863

UNIPROT ID:

Q86UL8

Full Name:

membrane associated guanylate kinase, WW and PDZ domain containing 2

Calculated MW:

1455 aa, 159 kDa

Purification Method:

Protein A purification

CloneNo.:

240512E9

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00529-1, MAGI2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83561-3-PBS. Detection antibody: 83561-4-PBS. Standard: Ag18279. Range: 0.781-100 ng/mL.