

For Research Use Only

# SYTL4 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83564-4-PBS

## Basic Information

Catalog Number: 83564-4-PBS	GenBank Accession Number: BC014913	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 94121	CloneNo.: 240556F9
Source: Rabbit	UNIPROT ID: Q96C24	
Isotype: IgG	Full Name: synaptotagmin-like 4	
Immunogen Catalog Number: AG2772	Calculated MW: 671 aa, 76 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00538-1, SYTL4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83564-3-PBS. Detection antibody: 83564-4-PBS. Standard: Ag2772. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00538-3, SYTL4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83564-4-PBS. Detection antibody: 83564-1-PBS. Standard: Ag2772. Range: 1.25-80 ng/mL