

For Research Use Only

HIBCH Recombinant antibody, PBS Only (Detector)



Catalog Number: 83569-4-PBS

Basic Information

Catalog Number: 83569-4-PBS	GenBank Accession Number: BC067822	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 26275	CloneNo.: 240590D9
Source: Rabbit	UNIPROT ID: Q6NVY1	
Isotype: IgG	Full Name: 3-hydroxyisobutyryl-Coenzyme A hydrolase	
Immunogen Catalog Number: AG6149	Calculated MW: 43 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00560-1, HibCH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83569-2-PBS. Detection antibody: 83569-4-PBS. Standard: Ag6149. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00560-2, HibCH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83569-3-PBS. Detection antibody: 83569-4-PBS. Standard: Ag6149. Range: 0.625-80 ng/mL