

For Research Use Only

ANAPC10 Recombinant antibody, PBS Only (Detector)



Catalog Number:83570-3-PBS

Basic Information

Catalog Number: 83570-3-PBS	GenBank Accession Number: BC005217	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10393	CloneNo.: 240492D8
Source: Rabbit	UNIPROT ID: Q9UM13	
Isotype: IgG	Full Name: anaphase promoting complex subunit 10	
Immunogen Catalog Number: AG8487	Calculated MW: 185 aa, 21 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00563-1, ANAPC10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83570-2-PBS. Detection antibody: 83570-3-PBS. Standard: Ag8487. Range: 0.313-40 ng/mL.